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warning

Please note the following points seriously when using the device and maintaining it:

- 1- Do not look directly into the LEDs with naked eyes. Light will cause temporary blindness.
- 2- Keep out of reach of children.
- 3- Refrain from putting too much pressure or hitting the device.
- 4- When cleaning And disinfecting, use a soft dampened cloth with disinfectant liquid.
- 5- Refrain from dropping any liquid on the device.
- 6- When the device is charging, do not use it
- 7- Do not leave the device in direct sunlight.
- 8- Do not leave the device in warm places.

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- 1- The less light, the better the functioning of the device. In case the vessel is not clearly visible, make the room darker.
- 2- In order to increase the lifetime of the device, highlight state should be used in cases where it is necessary.
- 3- In case you use high light without stopping for a long time, the device automatically turns off. In this case, turn off the device for five minutes. Then you can use the device again.
- 4- In case the battery runs out of power, the device automatically turns off. In this case, charge the device and then use it.
- 5- When you do not use the device, it is better to put it in charging mode.

Important notes for better using of the device

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Introduction

Since tissues in babies are soft and the organ's veins are internal, it is more difficult to find their veins. One of the most common problems nurses and doctors have that is associated with infants and child is how to give them shots or inject IV cannula. VAP company, by gathering experienced personnel and investigating different products that is used in various parts of the world, and also using advice from several pediatricians provided this vein finder called VeinTech. VeinTech, by utilizing Transillumination technology, has a considerable effect in revealing many vessels that pediatricians need to treat children. Thus, by using this device, the problem of finding veins in children and hemorrhage due to ruptured blood vessels during injection or taking blood from the children is significantly solved which always caused harm and injury to children. The device consists of three red LEDs which according to the situation, one of three LEDs are used to find veins. Also in terms of shape, size and amount of lighting has been designed such that the user can make the best use according to their need.

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Introduction

1- After using the LEDs for 2:30 minutes without putting LEDs on off mode, LEDs will change their status to off mode automatically to control the thermal status of VeinTech. Therefore, to use VeinTech again, one should wait until the sensors allow user to turn on the LEDs (Note: This time would be 5 minutes approximately)

2- For charging the device, the 1st two-state switch should be on off mode.

(Please pay attention:

If the 1st two-state switch is on, the device will not be charged.)

Introducing different parts of VeinTech

4- LEDs

It is the place where the light that determines the vessel is came from and is composed of three LEDs.

5- Socket to connect the charger to the device

Here we can connect the charger to the device.

6- Determining the battery charge status

By hooking the charger to the device, this part turns on and is on while the device is being charged. So this part turns off when the battery is full.

Introducing different parts of VeinTech

1- LED switch

When the first two-state switch is not off, by pressing this switch one or three LEDs turn/turns on. This switch has two parts that above part produce low light and the below part produces high light.

2- First two-state switch

This switch has two states: first state turns on one LED, second state turns on three LEDs.

3- Second two-state switch

This switch has two states, first state is off, second is on. When the device is charging, this switch must be at the first state.